



# APEX ASEPTIC 1422A TYVEK DRAWSTRING CLOSURE AUTOCLAVE POUCH

## Typical Application & Characteristics for APEX ASEPTIC 1422A Tyvek® Drawstring Closure Autoclave Pouch:

- Designed for steam sterilization
- Tyvek® on both sides for maximum breathability during sterilization
- Eliminates need to cut material and the use of autoclave tape
- Fast and consistent for equipment preparation area
- Creates secondary level of protection during aseptic presentation of equipment
- Heat-sealed seams eliminate sewing threads and greatly reduce particulates
- Hydrophobic
- Resists tearing and punctures
- Superior microbial barrier
- Lint-free – Category 1 (Helmke Drum)
- Manufactured and packaged in ISO certified controlled cleanroom environment
- Certificates of Conformance sent with shipment



PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A Tyvek®
Seal Strength	N/mm	F-88	>.88
Basis Weight	g/m <sup>2</sup>	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M2/24 HR	ASTM 398	1700
Hydrostatic Head	cm H <sub>2</sub> O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm <sup>2</sup> )	108
Puncture Resistance	N	EN 863	10.8
Melting Point	C°		135

APEX ASEPTIC ItemCode	1422A Tyvek® Drawstring Closure Pouch with Steam Indicator
C-T2BG0101UD	1" x 1"
C-T2BG0105UD	1" x 5"
C-T2BG00106UD	1" x 6"
C-T2BG0109UD	1" x 9"
C-T2BG0111UD	1" x 11"
C-T2BG0117UD	1" x 17"
C-T2BG0120UD	1" x 20"
C-T2BG0124UD	1" x 24"
C-T2BG0202UD	2" x 2"
C-T2BG02-503UD	2" x 3"
C-T2BG0205UD	2" x 5"
C-T2BG0208UD	2" x 8"
C-T2BG0211UD	2" x 11"
C-T2BG0304UD	3" x 4"
C-T2BG0305UD	3" x 5"
C-T2BG0306UD	3" x 6"
C-T2BG0308UD	3" x 8"
C-T2BG0309UD	3" x 9"
C-T2BG03109UD	3" x 10"
C-T2BG0404HD	4" x 4"
C-T2BG0407HD	4" x 7"
C-T2BG0409HD	4" x 9"
C-T2BG0412HD	4" x 12"
C-T2BG0425HD	4" x 25"
C-T2BG0430HD	4" x 30"
C-T2BG0443HD	4" x 43"
C-T2BG0448HD	4" x 48"
C-T2BG0606HD	6" x 6"
C-T2BG0608HD	6" x 8"
C-T2BG0612HD	6" x 12"

C-T2BG0617HD	6" x 17"
C-T2BG0623HD	6" x 23"
C-T2BG0633HD	6" x 33"
C-T2BG0638HD	6" x 38"
C-T2BG0643HD	6" x 43"
C-T2BG0662HD	6" x 62"
C-T2BG0710HD	7" x 10"
C-T2BG0717HD	7" x 17"
C-T2BG0727HD	7" x 27"
C-T2BG0808HD	8" x 8"
C-T2BG0809HD	8" x 9"
C-T2BG0812HD	8" x 12"
C-T2BG0817HD	8" x 17"
C-T2BG0824HD	8" x 24"
C-T2BG0831HD	8" x 31"
C-T2BG0843HD	8" x 43"
C-T2BG0856HD	8" x 56"
C-T2BG0865HD	8" x 65"
C-T2BG0872HD	8" x 72"
C-T2BG09-512HD	9.5" x 12"
C-T2BG1015HD	10" x 15"
C-T2BG1020HD	10" x 20"
C-T2BG1024HD	10" x 24"
C-T2BG1027HD	10" x 27"
C-T2BG1030HD	10" x 30"
C-T2BG1037HD	10" x 37"
C-T2BG1043HD	10" x 43"
C-T2BG1057HD	10" x 57"
C-T2BG1062HD	10" x 62"
C-T2BG1212HD	12" x 12"
C-T2BG1216HD	12" x 16"
C-T2BG1218HD	12" x 18"
C-T2BG1220HD	12" x 20"
C-T2BG1224HD	12" x 24"
C-T2BG1227HD	12" x 27"

APEX ASEPTIC ItemCode	1422A Tyvek® Drawstring Closure Pouch with Steam Indicator
C-T2BG1236HD	12" x 36"
C-T2BG1248HD	12" x 48"
C-T2BG1416HD	14" x 16"
C-T2BG1418HD	14" x 18"
C-T2BG1420HD	14" x 20"
C-T2BG1422HD	14" x 22"
C-T2BG1426HD	14" x 26"
C-T2BG1431HD	14" x 31"
C-T2BG1519HD	15" x 19"
C-T2BG1526HD	15" x 26"
C-T2BG1530HD	15" x 30"
C-T2BG1618HD	16" x 18"
C-T2BG15-519-5HD	15.5" x 19.5"
C-T2BG1616HD	16" x 16"
C-T2BG1624HD	16" x 24"
C-T2BG1627HD	16" x 27"
C-T2BG1630HD	16" x 30"
C-T2BG1632HD	16" x 32"
C-T2BG1638HD	16" x 38"
C-T2BG1642HD	16" x 42"
C-T2BG1544HD	15" x 44"
C-T2BG1819HD	18" x 19"
C-T2BG1824HD	18" x 24"
C-T2BG1830HD	18" x 30"
C-T2BG1836HD	18" x 36"
C-T2BG1848HD	18" x 48"
C-T2BG1854HD	18" x 54"
C-T2BG1855HD	18" x 55"
C-T2BG1870HD	18" x 70"
C-T2BG2020HD	20" x 20"
C-T2BG2022HD	20" x 22"

C-T2BG2024HD	20" x 24"
C-T2BG2026HD	20" x 26"
C-T2BG2030HD	20" x 30"
C-T2BG2033HD	20" x 33"
C-T2BG2037HD	20" x 37"
C-T2BG2040HD	20" x 40"
C-T2BG2218HD	22" x 18" 22" Side Open
C-T2BG2222HD	22" x 22"
C-T2BG2225HD	22" x 25"
C-T2BG2230HD	22" x 30"
C-T2BG2325HD	23" x 25"
C-T2BG2424HD	24" x 24"
C-T2BG2430HD	24" x 30"
C-T2BG2436HD	24" x 36"
C-T2BG2440HD	24" x 40"
C-T2BG2442HD	24" x 42"
C-T2BG2448HD	24" x 48"
C-T2BG2455HD	24" x 55"
C-T2BG2462HD	24" x 62"
C-T2BG2632HD	26" x 32"
C-T2BG2634HD	26" x 34"
C-T2BG2642HD	26" x 42"
C-T2BG2818HD	28" x 18" 28" Side Open
C-T2BG2827HD	28" x 27"
C-T2BG2834HD	28" x 34"
C-T2BG2838HD	28" x 38"
C-T2BG2844HD	28" x 44"
C-T2BG3030HD	30" x 30"
C-T2BG3036HD	30" x 36"
C-T2BG3040HD	30" x 40"
6-T2T2-3045HD	30" x 45"
6-T2T2-3053HD	30" x 53"

6-T2T2-3060HD	30" x 60"
C-T2BG3246HD	32" x 46"
C-T2BG3253HD	32" x 53"
C-T2BG3258HD	32" x 58"
C-T2BG3263HD	32" x 63"
C-T2BG3434HD	34" x 34"
C-T2BG3437HD	34" x 37"
C-T2BG3441HD	34" x 41"
C-T2BG3628HD	36" x 28"
C-T2BG3636HD	36" x 36"
C-T2BG3638HD	36" x 38"
C-T2BG3643HD	36" x 43"
C-T2BG3646HD	36" x 46"
C-T2BG3648HD	36" x 48"
C-T2BG3664HD	36" x 64"
C-T2BG4040HD	40" x 40"
C-T2BG4047HD	40" x 47"
C-T2BG4049HD	40" x 49"
C-T2BG4052HD	40" x 52"
C-T2BG4545HD	45" x 45"
C-T2BG4848HD	48" x 48"
C-T2BG4850HD	48" x 50"
C-T2BG4860HD	48" x 60"
C-T2BG4864HD	48" x 64"



**APEX ASEPTIC** | Pharmaceutical Cleanroom Division of Valumax  
 120 E. Gloucester Pike | Barrington, NJ 08007 | 856-805-7983  
[www.ApexAseptic.com](http://www.ApexAseptic.com)