



APEX ASEPTIC 1422A TYVEK® EQUIPMENT COVERS

ABOUT

Tyvek® Covers are designed to keep cleaned pharmaceutical equipment clean and protected.

FEATURES

- Individually poly bagged and labeled for easy identification
- Breathable Tyvek®
- Can incorporate elastic, zippers or snaps to ensure form fit and function
- Allows residual moisture from cleaning process to evaporate
- Hydrophobic
- Resists tearing and punctures
- Anti-static
- Lint-free – Category 1 (Helmke Drum)
- Manufactured and packaged in ISO certified controlled cleanroom environment
- Certificates of Conformance sent with shipment



PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A Tyvek®
Basis Weight	g/m ²	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M2/24 HR	ASTM 398	1700
Hydrostatic Head	cm H ₂ O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm ²)	108
Puncture Resistance	N	EN 863	10.8
Melting Point	C°		135

Ordering Information for Stocked Sizes:

APEX ASEPTIC ItemCode	Tyvek® Cover
<i>Rectangular Equipment Covers:</i>	
C-T2CL242424	24" x 24" x 24"
C-T2CL322412	32" x 24" x 12"
C-T2CL322824	32" x 28" x 24"
C-T2CL362442	36" x 24" x 42"
C-T2CL363660	36" x 36" x 60"
C-T2CL423460	42" x 34" x 60"
C-T2CL443532	44" x 35" x 32"
C-T2CL443532	44" x 44" x 60"
C-T2CL484860	48" x 48" x 60"
C-T2CL482430	48" x 24" x 30"
C-T2CL505014	50" x 50" x 14"
C-T2CL505050	50" x 50" x 50"
C-T2CL561444	56" x 14" x 44"
C-T2CL603072	60" x 30" x 72"
C-T2CL603460	60" x 34" x 60"
C-T2CL604444	60" x 44" x 44"
C-T2CL604848	60" x 48" x 48"
C-T2CL605060	60" x 50" x 60"
C-T2CL607860	60" x 78" x 60"

APEX ASEPTIC ItemCode	Tyvek® Cover
<i>Cylindrical Equipment Covers:</i>	
C-T2CL50D50	50" x 50
C-T2CL50D60	50" x 60
C-T2CL60D60	60" x 60
C-T2CL14D60	14" x 60
C-T2CL20D48	20" x 48"
C-T2CL32D60	32" x 60
C-T2CL48D36	48" x 36"
C-T2CL40D32	40" x 32"
C-T2CL72D72	72" x 72"

***All Equipment Covers can be manufactured with elastic or zippers.**

To order covers with elastic:
Add "E" to the end of any of the item codes.

To order covers with zippers:
Add "Z" to the end of any of the item codes.



APEX ASEPTIC | Pharmaceutical Cleanroom Division of Valumax
120 E. Gloucester Pike | Barrington, NJ 08007 | 856-805-7983
www.ApexAseptic.com