



APEX ASEPTIC 1422A TYVEK® STOPPER BOWL COVERS

ABOUT

Used as an autoclavable cover these Apex Aseptic 1422A Stopper Bowl Covers eliminate the need for time consuming wrapping with sheets of material and autoclave tape. When used in conjunction with a Tyvek® Bag it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover remains in place and provides protection from microbial and particulate contamination while remaining set-up activities are performed.

FEATURES

- Steam Sterilization Indicator
- Lint-free – Category 1 (Helmke Drum)
- Specialized tunneled sewing technique eliminates exposure of elastic and threads - greatly reducing particulate contamination risk
- Autoclavable at 121C-125C
- Hydrophobic
- Resists tearing and punctures
- Manufactured and packaged in ISO certified controlled cleanroom environment
- Conveniently packaged: 2 per heat-sealed labeled easy-tear poly bag
- Certificates of Conformance sent with shipment

STEAM STERILIZATION INDICATOR

SPECIALIZED TUNNELED SEWING

REDUCED COST & DELIVERY TIME

BUY DIRECT FROM MANUFACTURER

AVAILABLE FOR IMMEDIATE SHIPMENT
FROM LOCAL WAREHOUSES



PROPERTY	UNITS	TEST METHOD	TYPICAL VALUES 1422A Tyvek®
Seal Strength	N/mm	F-88	>.88
Basis Weight	g/m ²	EN ISO 536	41
Total Thickness	μ	EN 20534	130
Air Permeability	s	ISO 5636-5	20
Water Vapor Permeability	G/M2/24 HR	ASTM 398	1700
Hydrostatic Head	cm H ₂ O	ISO 811	120
Trapezoidal Tear MD/CD	N	ISO 9073-4	26.1/30.6
Linting Propensity	Rating	BS 6909 (Shirley Method 21)	Pass Rating: Excellent
Burst Resistance	kPa	ISO 2960 (50 cm ²)	108
Puncture Resistance	N	EN 863	10.8
Melting Point	C°		135

Ordering Information for Stocked Sizes

APEX ASEPTIC ItemCode	1422A Tyvek® Bowl Cover With Elastic and Steam Indicator
C-T2R20AST2	20", 2 per Sealed Poly Bag
C-T2R24AST2	24", 2 Per Sealed Poly Bag
C-T2R30AST2	30", 2 Per Sealed Poly Bag
C-T2R36AST2	36", 2 Per Sealed Poly Bag
C-T2R42AST2	42", 2 per Sealed Poly Bag
C-T2R48AST2	48", 2 Per Sealed Poly Bag
C-T2R60AST2	60", 2 Per Sealed Poly Bag

CONTACT

APEX ASEPTIC

Pharmaceutical Cleanroom Division of Valumax

120 E. Gloucester Pike
Barrington, NJ 08007
856-805-7983

www.ApexAseptic.com

