



APEX ASEPTIC NON-WOVEN STOPPER BOWL COVERS

ABOUT

These Apex Aseptic Non-Woven Stopper Bowl Covers eliminate the need for time consuming wrapping with sheets of material and autoclave tape. When used in conjunction with an autoclavable bag it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover remains in place and provides protection from microbial and particulate contamination while remaining set-up activities are performed.

REDUCED COST & DELIVERY TIME

BUY DIRECT FROM MANUFACTURER

AVAILABLE FOR IMMEDIATE SHIPMENT
FROM LOCAL WAREHOUSES

FEATURES

- Specialized tunneled sewing technique eliminates exposure of elastic and threads - greatly reducing particulate contamination risk
- Steam Sterilization Indicator
- Manufactured and packaged in ISO certified controlled cleanroom environment
- Autoclavable up to 135C
- Certificates of Conformance sent with shipment
- Conveniently packaged: 2 per heat-sealed, labeled easy-tear poly bag
- Colors: Blue



PROPERTY	UNITS	NON-WOVEN WRAPPING MATERIAL
Basis Weight	Gm/m ²	61
Dry Tensile Strength	N/m	1790
Elongation Dry	%	27.3
Bursting Strength	kPa	300

Ordering Information for Stocked Sizes

APEX ASEPTIC ItemCode	BLUE Non-woven Stopper Bowl Cover With Steam Indicator
C-R20NW- B	20"
C-R24NW- B	24"
C-R30NW- B	30"
C-R36NW- B	36"
C-R42NW- B	42"
C-R48NW- B	48"
C-R52NW- B	52"
C-R60NW- B	60"



CONTACT

APEX ASEPTIC

Pharmaceutical Cleanroom Division of Valumax

120 E. Gloucester Pike
 Barrington, NJ 08007
 856-805-7983

www.ApexAseptic.com