



APEX ASEPTIC MEDICAL GRADE PAPER STOPPER BOWL COVERS

ABOUT

These Apex Aseptic Stopper Bowl Covers eliminate the need for time consuming wrapping with sheets of material and autoclave tape. When used in conjunction with an autoclavable bag it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover remains in place and provides protection from microbial and particulate contamination while remaining set-up activities are performed.

REDUCED COST & DELIVERY TIME

BUY DIRECT FROM MANUFACTURER

AVAILABLE FOR IMMEDIATE SHIPMENT
FROM LOCAL WAREHOUSES

FEATURES

- Specialized tunneled sewing technique eliminates exposure of elastic and threads - greatly reducing particulate contamination risk
- Steam Sterilization Indicator
- Manufactured and packaged in ISO certified controlled cleanroom environment
- Autoclavable up to 135C
- Certificates of conformance sent with shipment
- Conveniently packaged: 2 per heat-sealed, labeled easy-tear poly bag
- Colors: Blue, Green



PROPERTY	UNITS	BLUE & GREEN MEDICAL GRADE CREPE PAPER
Basis Weight	Gm/m ²	60 (+/-5)
Dry Tensile Strength MD	N	>15
Dry Tensile Strength CD	N	>10
Wet Tensile Strength MD	N	>5
Elongation MD	%	>5
Elongation CD	%	>2

Ordering Information for Stocked Sizes

APEX ASEPTIC ItemCode	BLUE Medical Grade Paper Stopper Bowl Cover w/ Steam Indicator
C-R20CP- B	20"
C-R24CP- B	24"
C-R30CP- B	30"
C-R36CP- B	36"
C-R42CP- B	42"
C-R48CP- B	48"
C-R52CP- B	52"
C-R60CP- B	60"

APEX ASEPTIC ItemCode	GREEN Medical Grade Paper Stopper Bowl Cover w/ Steam Indicator
C-R20CP- G	20"
C-R24CP- G	24"
C-R30CP- G	30"
C-R36CP- G	36"
C-R42CP- G	42"
C-R48CP- G	48"
C-R52CP- G	52"
C-R60CP- G	60"

CONTACT

APEX ASEPTIC

Pharmaceutical Cleanroom Division of Valumax

120 E. Gloucester Pike
Barrington, NJ 08007
856-805-7983

www.ApexAseptic.com

