



APEX ASEPTIC BLUE NONWOVEN TINY MUFFS

REDUCED COST & DELIVERY TIME

BUY DIRECT FROM MANUFACTURER

AVAILABLE FOR IMMEDIATE SHIPMENT
FROM LOCAL WAREHOUSES

ABOUT

Used as an autoclavable cover these Apex Aseptic Nonwoven Tiny Muffs eliminate the need for time consuming wrapping with sheets of material and autoclave tape.

FEATURES

- Manufactured and packaged in ISO certified controlled cleanroom environment
- Autoclavable up to 135C
- Certificates of Conformance sent with shipment
- Conveniently packaged: 100 per heat-sealed, labeled easy-tear poly bag
- Colors: Blue



PROPERTY	UNITS	Nonwoven Wrapping Material
Basis Weight	Gm/m ²	61
Dry Tensile Strength	N/m	1790
Elongation Dry	%	27.3
Bursting Strength	kPa	300

Ordering Information for Stocked Sizes

APEX ASEPTIC ItemCode	BLUE Nonwoven Tiny Muffs
C-R02NW- B	2" used for openings ½" and smaller
C-R03NW- B	3" used for openings ½" to 1"
C-R04NW- B	4" used for openings 1" to 1 ¾"
C-R05NW- B	5" used for openings 1 ¾" to 2 ½"
C-R07NW- B	7" used for openings 2 ½" to 4"
C-R08NW- B	8" used for openings 3 ½" to 4 ¾"
C-R09NW- B	9" used for openings 4" to 5 ½"
C-R012NW- B	12" used for openings 5 ½" to 8"
C-R016NW- B	16" used for openings 8" to 12"
C-R018NW- B	18" used for openings 10" to 14"



CONTACT

APEX ASEPTIC

Pharmaceutical Cleanroom Division of Valumax

120 E. Gloucester Pike
Barrington, NJ 08007
856-805-7983

www.ApexAseptic.com