



# APEX ASEPTIC MEDICAL GRADE PAPER TINY MUFFS

REDUCED COST & DELIVERY TIME

BUY DIRECT FROM MANUFACTURER

AVAILABLE FOR IMMEDIATE SHIPMENT  
FROM LOCAL WAREHOUSES

## ABOUT

Used as an autoclavable cover these Apex Aseptic Medical Grade Paper Tiny Muffs eliminate the need for time consuming wrapping with sheets of material and autoclave tape.

## FEATURES

- Manufactured and packaged in ISO certified controlled cleanroom environment
- Autoclavable up to 135C
- Certificates of Conformance sent with shipment
- Conveniently packaged: 100 per heat-sealed, labeled easy-tear poly bag
- Colors: Blue, Green



PROPERTY	Units	Green & Blue Medical Grade Crepe Paper
Basis Weight	Gm/m <sup>2</sup>	60 (+/-5)
Dry Tensile Strength MD	N	>15
Dry Tensile Strength CD	N	>10
Wet Tensile Strength MD	N	>5
Elongation MD	%	>5
Elongation CD	%	>2

## Ordering Information for Stocked Sizes

APEX ASEPTIC ItemCode	BLUE Medical Grade Paper Tiny Muffs
C-R02CP-B	2" used for openings ½" and smaller
C-R03CP-B	3" used for openings ½" to 1"
C-R04CP-B	4" used for openings 1" to 1 ¾"
<i>2" - 4" (100 per sealed poly bag, 30 poly bags per case)</i>	
C-R05CP-B	5" used for openings 1 ¾" to 2 ½"
<i>5" (100 per sealed poly bag, 20 poly bags per case)</i>	
C-R07CP-B	7" used for openings 2 ½" to 4"
C-R08CP-B	8" used for openings 3 ½" to 4 ¾"
C-R09CP-B	9" used for openings 4" to 5 ½"
C-R12CP-B	12" used for openings 5 ½" to 8"
<i>7" - 12" (100 per sealed poly bag, 10 poly bags per case)</i>	
C-R16CP-B	16" used for openings 8" to 12"
C-R18CP-B	18" used for openings 10" to 14"
<i>16" - 18" (50 per sealed poly bag, 10 poly bags per case)</i>	

APEX ASEPTIC ItemCode	GREEN Medical Grade Paper Tiny Muffs
C-R02CP-G	2" used for openings ½" and smaller
C-R03CP-G	3" used for openings ½" to 1"
C-R04CP-G	4" used for openings 1" to 1 ¾"
<i>2" - 4" (100 per sealed poly bag, 30 poly bags per case)</i>	
C-R05CP-G	5" used for openings 1 ¾" to 2 ½"
<i>5" (100 per sealed poly bag, 20 poly bags per case)</i>	
C-R07CP-G	7" used for openings 2 ½" to 4"
C-R08CP-G	8" used for openings 3 ½" to 4 ¾"
C-R09CP-G	9" used for openings 4" to 5 ½"
C-R12CP-G	12" used for openings 5 ½" to 8"
<i>7" - 12" (100 per sealed poly bag, 10 poly bags per case)</i>	
C-R16CP-G	16" used for openings 8" to 12"
C-R18CP-G	18" used for openings 10" to 14"
<i>16" - 18" (50 per sealed poly bag, 10 poly bags per case)</i>	

## CONTACT

### APEX ASEPTIC

Pharmaceutical Cleanroom Division of Valumax

120 E. Gloucester Pike  
Barrington, NJ 08007  
856-805-7983

[www.ApexAseptic.com](http://www.ApexAseptic.com)

